



922290@swan.ac.uk

Investigating charge carrier recombination in perovskite materials

Emmanuel V. Péan, Stoichko Dimitrov, Catherine S. De Castro and Matthew L. Davies

Next Generation Materials for Solar Photovoltaics 2020

Thanks to their ideal band-gap and high absorption coefficient, perovskite materials are good candidates for the next generation of solar cells with an impressive certified 25.2% power conversion efficiency¹. In particular, $\text{CH}_3\text{NH}_3\text{PbI}_3$ (MAPI) solar cells presents a simple solution processing alongside an easy deposition method by spin coating. However, point defects such as iodide vacancies have been shown to be responsible for degradation and performance losses in such devices [2-3]. Time resolved photoluminescence (TRPL) allows to investigate charge carrier recombinations and defects in perovskites. Here the interpretation of such measurement is discussed.

- ▶ Defect-induced trap states are responsible for degradation and performance losses in perovskites [2-3];
- ▶ TRPL is a measure of the probability of a photon to be emitted by the perovskite through a radiative charge carrier recombination after a certain time after a short excitation pulse.
- ▶ Radiative recombinations compete with other recombination processes including trap states mediated recombinations.
- ▶ Time correlated single photon counting (TCPSC) allows to measure TRPL and thus to investigate trap states in perovskites.




Figure 1: Insides of a TCSPC. A laser repeatedly excites electrons in the sample and after each excitation pulse, 0 or 1 emitted photons are measured.





Figure 2: Example of photons measured after 5 consecutive pulses at $t=0$ (top) and TRPL intensity obtained after multiple excitation pulses (bottom).

Introduction


n_e : Photoexcited electron concentration
 n_h : Photoexcited hole concentration
 n_t : Trapped electron concentration
 N_T : Deep trap state concentration




Figure 3: Evolution of the charge carrier concentrations and TRPL intensity after 3 consecutive excitation pulses for excitation pulse repetition rates of a) 20 μs and b) 5 μs simulated using constants from [5].




Figure 4: a) Fitting of TRPL decays of a MAPI thin film on glass measured at different excitation fluencies and b) results of the fitting of the same data using a wide range of guess parameters.

Conclusion

- ▶ Charge carrier recombination in perovskites require two carriers to happen and thus depend on the carrier concentrations;
- ▶ It is possible to fit experimental TRPL decays by simulating a single excitation pulse;
- ▶ However, this requires to ascertain that all experimental pulses lead to the same variation of the carrier concentrations;
- ▶ Samples need to be measured with a wide range of excitation fluencies in order to observe all the recombination processes.

References

- [1] National Renewable Energy Laboratory, "Best Research-Cell Efficiency Chart, NREL, 2019.
 [2] R. Brenes et al, *Adv. Mater.* 2018, 30, 1706208.

- [3] D. W. DeQuilettes et al, *Nat. Commun.* 2016, 7, 11683.
 [4] E. M. Hutter et al, *J. Phys. Chem. Lett.* 2015, 6, 3082.
 [5] R. Brenes et al, *Joule* 2017, 1, 155.

Acknowledgments

This work has been funded by the Zienkiewicz Scholarship, the Sêr Cymru Solar initiative, the EPSRC (EP/R016666/1, EP/S001336/1 and EP/M028267/1) and the European Regional Development Fund through the Welsh Government

Led by

 Swansea University
Prifysgol Abertawe

Funders
